ReStore™ Microcatheter

Manufacturer	European Authorized Representative
Reverse Medical Corp.	MediMark Europe Sarl
13900 Alton Parkway	11, rue Emile Zola. BP 2332
Suite 123	F-38033 Grenoble Cedex 2
Irvine, CA 92618	France

Figure 1. ReStore™ Microcatheter



Figure 2. ReAct™ Stylet

Device Description

The ReStore™ Microcatheter consists of a flexible, tapered reinforced composite catheter with a braided mesh retrieval element attached to the distal region of the catheter. The inner lumen can accommodate a guidewire to aid in placement of the catheter. The catheter contains radiopaque markers to facilitate fluoroscopic visualization— a distal tip marker, markers at the proximal and distal ends of the retrieval element, and circumferentially placed markers within the retrieval element. The retrieval element is deployed through the advancement of a guidewire or the ReAct™ stylet through the lumen to radially expand the retrieval element. The selection of guidewire sizes in accordance with device compatibility allows for guidewires to be used for either retriever placement or for deployment of the element.

The proximal end of the catheter has a luer fitting to allow attachment of accessories and infusion of liquids through the catheter. The catheter is coated with a hydrophilic coating. The catheter is designed to be used with a guide catheter.

The catheter is offered in various retrieval element sizes to accommodate physician preferences.

The ReStore™ Microcatheter can be used in vessel diameters ranging from 2mm to 4mm.

Each ReStore™ Microcatheter may be packaged with a ReAct™ stylet and stylet handle. The ReAct™ stylet is a stainless steel wire with a platinum tip that may be used to activate the ReStore™ Microcatheter's retrieval element. The ReAct™ stylet is provided in the same sterile pouch as the ReStore™ Microcatheter.

Indications for Use

The ReStore™ Microcatheter is intended to restore blood flow in the neurovasculature by removing thrombus in patients experiencing ischemic stroke. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV-t-PA therapy are candidates for treatment with ReStore™

microcatheters are also indicated for use in the retrieval of foreign bodies misplaced during interventional radiological procedures in the neuro, peripheral and coronary vasculature.

Potential Complications

Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with possible complications which may occur during or after the procedure. Possible complications include, but are not limited to the following: air embolism; hematoma or hemorrhage at puncture site; infection; distal embolization; vessel spasm, thrombosis, dissection, or perforation; emboli; acute occlusion; ischemia; intracranial hemorrhage; false aneurysm formation; neurological deficits including stroke; and death.

Precaution: To reduce the probability of coating damage in tortuous vasculature, wide lumen guide catheters are recommended for use with the ReStore™ Microcatheter.

Compatibility

The ReStore™ Microcatheter is compatible with the following catheters and guidewires:

- Guide catheters with minimum internal diameter 0.052 in.
- Neuro guidewires for navigation and activation as specified in table below.

Compatibility and product specifications:									
Model	Navigation Guidewire	Activation Guidewire	A (cm)	B (mm)	C (mm)	D (mm)	E (mm)	F (mm)	G (cm)
RST114-3.0x15	0.014 in.	0.016 in.	160	3.0	8	15	1	10	180
RST114-3.5x20	0.014 in.	0.016 in.	160	3.5	8	20	1	15	180
RST114-4.0x25	0.014 in.	0.016 in.	160	4.0	8	25	1	20	180

Recommended Preparation Procedure

Continuous Flush

Precaution: In order to achieve optimal performance of the ReStore™ Microcatheter, it is critical that a continuous infusion of appropriate flush solution be maintained between a) the guide catheter and the ReStore™ Microcatheter, and b) the ReStore™ Microcatheter and the guidewire during a procedure. Flushing aids in preventing contrast crystal formation and/or clotting on both the intraluminal device and inside the catheter lumen.

The recommended continuous flush setup as illustrated in Figure 3, requires two stopcocks, two rotating hemostatic valves (RHV) (Tuohy-Borst type), two flush lines and two 3 cc syringes. The RHVs provide a fluid tight seal and are attached to the guide catheter and ReStore™ Microcatheter hub. The stopcocks attach to the RHV sidearm which becomes infusion ports for appropriate flush or contrast media injection.

ReStore™ Microcatheter

Flush Set up to ReStore™ Microcatheter

1. Connect a drip line for the appropriate flush solution between the ReStore™ Microcatheter and the guidewire using the stopcock attached to the 1st RHV on the on the ReStore™ Microcatheter. In addition to providing a flushing port, the RHV provides a fluid tight seal between the ReStore™ Microcatheter and guidewire or ReAct™ stylet.
2. Prepare the ReStore™ Microcatheter for use by flushing with appropriate flush solution using a 3cc syringe.
3. Connect a drip line for the appropriate flush solution between the guide catheter and the ReStore™ Microcatheter using the stopcock attached to the 2nd RHV on the on the guide catheter. In addition to providing a flushing port, the RHV provides a fluid tight seal between the guide catheter and the ReStore™ Microcatheter.
Precaution: In order to achieve optimal performance of the ReStore™ Microcatheter, it is critical that a continuous infusion of appropriate flush solution be maintained between a) the guide catheter and the ReStore™ Microcatheter, and b) the ReStore™ Microcatheter and guidewires or stylet during the procedure. Flushing aids in preventing contrast crystal formation and/or clotting on the guidewire or stylet and inside the catheter lumen. Inadequate flow can cause serious injury to the patient.

Warning

The ReStore™ Microcatheter is compatible for use with contrast media, Heparin and Normal saline (0.9% Sodium Chloride Injection) at room temperature only. Use of other chemicals is not recommended.

Precaution: Check that all fittings are secure so that air is not introduced into a guide catheter or infusion during continuous flush.

Directions for Use

1. Place the appropriate guide catheter according to the guide catheter manufacturer's instructions for use.
2. Carefully insert the ReStore™ Microcatheter into 2nd RHV connected to the guide catheter and advance the ReStore™ Microcatheter through the RHV of the guide catheter. Take care not to kink the ReStore™ Microcatheter shaft. Tighten the RHV o-ring valve just enough to prevent back flow, but not enough to hinder movement of the ReStore™ Microcatheter.
Precaution: Free movement of the guidewire or ReAct™ stylet within the ReStore™ Microcatheter is an important feature of the system. Whenever possible, test the system for resistance prior to use, and replace device if necessary.

Warning

Never advance, rotate or withdraw an endovascular device against resistance until cause of resistance is determined. Movement or torquing of ReStore™ Microcatheter or guidewire or ReAct™ stylet against resistance could dislodge a clot, perforate a vessel wall, cause vasospasm, or damage ReStore™ Microcatheter or guidewire or ReAct™ stylet. In severe cases, tip separation of ReStore™ Microcatheter or guidewire may occur and this can cause vascular damage or other patient injury.

3. Using conventional catheterization techniques and fluoroscopic guidance, place the ReStore™ Microcatheter into the selected target vessel by carefully manipulating the navigation guidewire within the ReStore™ Microcatheter.

4. For retrieval of a foreign body position the ReStore™ Microcatheter with the retrieval element directly adjacent to and as close as possible to the object or foreign body to be retrieved. For retrieval of a thrombus or clot position the ReStore™ Microcatheter with the retrieval element distal to and as close as possible to the thrombus or clot to be retrieved.
5. Once positioned, remove the navigation guidewire from the ReStore™ Microcatheter.
6. Insert the activation guidewire or ReAct™ stylet into the ReStore™ Microcatheter until the stylet handle approaches and is 1cm from the 1st RHV.
7. Slowly advance the stylet handle until the positioning marker bands are visualized on fluoroscopy and their movement towards each other is observed. Radial expansion of the element must be confirmed fluoroscopically by observing the radial expansion of four marker bands positioned at the center of the element. The element diameter should remain slightly smaller than the vessel diameter.
8. For foreign body retrieval the target object should be captured within the expanded element by SLOWLY withdrawing the activation guidewire or ReAct™ stylet handle to close the retrieval element and entrap the foreign object. For thrombus retrieval the retrieval element should remain in the expanded state.
9. With foreign body or thrombus securely captured SLOWLY and carefully withdraw the ReStore™ Microcatheter and ReAct™ stylet together with the object or clot captured to the guide catheter tip.
10. Apply vigorous aspiration to the guide catheter using 60 ml syringe and withdraw the ReStore™ Microcatheter and the object or clot captured into the guide catheter and withdraw the whole system including the guide catheter.
Precaution: Do not attempt to remove the ReStore™ Microcatheter with snared object by pulling through the guide catheter. The object may catch on the end of the guide catheter, separate from the retrieval element and return to the vasculature.
11. Inspect the ReStore™ Microcatheter for damage to the element or catheter body. If not damaged, the ReStore™ Microcatheter may be used for up to two (2) retrieval attempts. A retrieval attempt is one advancement and complete withdrawal cycle.
Precaution: To reduce the risk of vessel damage, adhere to the following recommendations:
- Take care to expand the ReStore™ Microcatheter to the appropriate vessel diameter at the intended site of deployment- Do not perform more than six (6) retrieval attempts in same vessel using the ReStore™ Microcatheters.

ADDITIONAL WARNINGS AND PRECAUTIONS

Precaution: Observe all element movement in the vessels using fluoroscopy and contrast media.

Warning

Before retrieval element is opened or moved, guidewire or ReAct™ stylet movement should be carefully monitored under fluoroscopy to prevent the possibility of vessel damage.

Warning

Never advance or withdraw an intraluminal device against resistance until the cause of resistance is determined by fluoroscopy.

ENGLISH

INSTRUCTIONS FOR USE

ReStore™ Microcatheter

SA03013 Rev. P 4 of 14

Warning

If flow through the ReStore™ Microcatheter becomes restricted, do not attempt to clear catheter lumen by infusion; remedy the cause of the blockage, or replace the blocked ReStore™ Microcatheter before resuming procedure.

Warning

Pressure in excess of 100 psi may result in ruptured microcatheter.


Warning

Dispose of used devices in accordance with established protocols for handling bio-hazardous materials.

PACKAGING & STORAGE

- The ReStore™ Microcatheter and ReAct™ stylet are sterile and non-pyrogenic in unopened, undamaged packaging.
- Do not use open or damaged packages.
- Do not use if any portion of the ReStore™ Microcatheter or ReAct™ stylet appears damaged.
- The Reverse Medical ReStore™ Microcatheter is for single use only and should not be cleaned re-processed, resterilized or reused as these processes may damage the components and surface of the device and compromise performance. In addition, the risks of infection of a reprocessed device have not been qualified to reduce this risk to an acceptable level. Structural integrity and/or function may be impaired by reuse or cleaning.
- The ReStore™ Microcatheter and the ReAct™ stylet are delicate instruments and should be handled carefully. Do not use a device that shows signs of damage. Damage may prevent device from functioning and may cause complications.
- Use by "Use By" date.
- Exposure to temperatures above 54°C (130°F) may damage device and accessories. Do not autoclave.
- Do not expose the devices to solvents.

Symbols Glossary

i	Consult instructions for use	
Y	Caution, consult accompanying documents	
D	Do not reuse	B Do not re-sterilize
		h Catalogue Number
L	Do not use if package is damaged	
IQ	Sterilized using Ethylene Oxide	 Use By
W	Keep away from sunlight	P Authorized Representative in the European Community
 Batch Code		 Contents
 Nonpyrogenic		M Manufacturer
 Rx Only	Prescription only- device restricted to use by or on order of a physician	
	Flush package hoop with saline before removing device	
CE0297		S 54°C Upper limit of temperature 54°C

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